

Congratulations on your purchase of TiFinity Toothbrush.

Proper brushing promotes healthy gums that will literally save thousands of dollars in future dental bills. Allow me to offer some brushing tips:

- Gently tease the tips of the bristles towards the gum line so that you concentrate on cleaning at the neck of the teeth and below the gums. Use small circular motion with your TiFinity Toothbrush.
- It is recommended you brush a minimum of twice daily for approximately two minutes each time. .
- Rinse your TiFinity Toothbrush off with hot water after every use.

As a practicing Dentist for over 28 years, most which was specialized in dental implants; It is my opinion you just purchased the most effective Home Care dental hygiene tool on the market. You will notice the amazing reach of the micro fine bristles; cleaning areas between the teeth and under the gums that simply cannot be reached with traditional toothbrushes. Once you get accustomed to TiFinity's super-elastic titanium bristles, you won't ever want to go back to nylon. It is our hope that TiFinity Toothbrush will raise your level of dental health to heights you never experience before and your teeth will stay healthy for the duration of your life.

Thank You for your Purchase  
Sincerely,  
Dr. Dane Robinson, President

[www.tifinitytoothbrush.com](http://www.tifinitytoothbrush.com)

